




BioXp™ Loading Map and Checklist– Cloning

Each BioXp™ Cloning kit includes Module A (+4°C), Module B (–20°C), and Module C (–20°C).

- Load reagents for a new run according to the following instructions:
 - Use appropriately sized reagents for your job; reagent sizes are denoted by stickers on boxes, e.g. **S** or **L**
 - If the door is closed, select "Unlock Door" from the instrument LCD screen and open the door
- Thaw –20°C components as directed below:
 - DNA Assembly Reagent Plate** (30 minutes at room temperature or 1 hour on ice)
 - GA Cloning Strip** and **Vector Strip** (30 minutes on ice)
- Load tips by aligning the tip tray notch with the **upper left corner** of each Tip Tray Retainer
 - Load 3 x 50 µL tips
 - Load 1 x 200 µL tips
- Add a minimum of 12 mL freshly prepared 70% ethanol to the reusable **Ethanol Reservoir**
 - Load **Ethanol Reservoir** in the right-most Reservoir Retainer position of the instrument deck

Note: Do not discard the **Ethanol Reservoir** after the run; keep for future use
- Load plates stored at 4°C:
 - Load the **Recovery Plate** onto the Recovery Chiller with the notch in the **upper left corner**
 - Load the **Oligo Vault™ Plate** into the Thermocycler with the notch in the **upper left corner**
- Briefly spin the strips. Load in the order listed and shown below, with the strip pinhole closest to you
 - Load the black **DNA Purification Strip** into position #1
 - Leave position #2 empty
 - Load the purple **GA Cloning Strip** into position #3
 - Load the red **pUCGA 1.0** or **Custom Vector Strip** into position #4
- Secure strips with spring-loaded arms while holding the strips in place.
- Spin the thawed **DNA Assembly Reagent Plate** for 1 minute at 500 rpm.
 - Load **DNA Assembly Reagent Plate** onto Reagent Chiller, notch in the **lower left corner**

Note: Be certain that the plate is properly seated within the chiller.
- Refer to the photo in the lower left panel below. Confirm that components are securely seated. Close the door.
- After the deck inspection ends, press **Start Now** or **Delay Start** to begin the run.



Guidelines to prepare the BioXp™ custom vector strip


1. Adjust the vector concentration according to the table:

Vector Size (kb)	Concentration (ng/μL)
3–5	15–20
5–7	20–25
7–9	25–30
9–12	30–35

2. Determine the volume of prepared linear vector you will add to wells A, C, E, and G of a BioXp™ Vector strip.


Number of BioXp™ Cloning Reactions	Volume
≤16 Reactions	12 μL
>16 Reactions	18 μL

3. Add 12 or 18 μL of prepared linear vector to wells A, C, E and G of a BioXp™ Vector strip. Do not seal the strip. Load the prepared custom strip on the deck with the pinhole closest to the front of the instrument..



EXAMPLE: To prepare the vector strip for 8 cloning reactions, prepare a 10 kb vector at a concentration of 30–35 ng/μL. Add 12 μL of prepared vector to the four strip wells (A, C, E and G). Total amount of required vector = 1.44 to 1.68 μg.

Final DNA product location

Deck at run completion	Items to keep		
	DNA product	Located within	Storage instructions
	Cloning Reaction (wells A1–H4)	Recovery Plate	Seal and store at 2°C to 8°C for up to one week or at –20°C for up to a year.
	Uncloned Tiles (wells A5–H8)	Oligo Vault™ Plate (foil covered)	Seal and store at 2°C to 8°C for up to one week.
	Ethanol Reservoir (empty and rinse for next use)		

Tips

After the BioXp™ run is complete, cloning reactions are ready for transformation. We recommend diluting cloning reactions 1:2 with molecular biology grade water before transformation.

Analyze Tiles– We recommend evaluating the success of the assembly reaction by running a gel containing the uncloned BioXp™ Tiles from the Oligo Vault™ Plate before transforming clones

Transformation– We recommend using *E. coli* 10G Chemically Competent Cells (Lucigen Cat. No. 60107, free samples available at Lucigen.com) or TransforMax™ EPI300™ Electrocompetent *E. coli* (Lucigen Cat. No. EC300110). If other cells are used, be certain to use competent cells with a transformation efficiency $\geq 1 \times 10^9$ CFU/μg pUC19.

Additional product information is available at www.sgidna.com/bioxp

Technical Services techservices@sgidna.com

Trademarks E. cloni® is a registered trademark of Lucigen Corporation.
TransforMax™ and EPI300™ are trademarks of Epicentre Technologies Corporation.
BioXp™, Oligo Vault™, and Bio360™ are trademarks of SGI-DNA.